OSTEO-X®

EN INSTRUCTIONS FOR USE

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SV BRUKSANVISNING

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TR KULLANIM KILAVUZU

EN INSTRUCTIONS FOR USE

DESCRIPTION

The Osteo-x is an adjustable external fixator which is used for the treatment of small bone fractures of the hand. The device is supplied in a fully assembled, sterile state; comprising of two K-wire blocks each containing two K-wires. The K-wire blocks are linked by a threaded shaft and two other support struts which provide structural integrity. The threaded shaft is designed such that the included adjustment tool can be used to alter the distance between the K-wire blocks. This can initially be used so that the device is sized for the patient and can be used subsequently to achieve distraction or compression.

GENERAL INSTRUCTIONS

The Osteo-x should be used for the treatment of fracture fixation and distraction of small hand bones.

When using the device for the distraction of bones, it is important to know that one full turn of the threaded shaft will result in a change of 0.5mm displacement between the two blocks.

INDICATIONS

The Osteo-x external fixator is designed for fracture fixation and distraction of bones in the hand/wrist.

- Small hand and wrist bone fractures
- Limb lengthening

CONTRAINDICATIONS

- Infection
- Insufficient bone quantity
- Poor bone quality/osteoporosis
- Allergic reaction to device materials
- Un-cooperative patient
- Potentially successful conservative treatment

The judgement of whether an implant is suitable for a particular patient is at the surgeon's discretion.

ADVERSE EFFECTS

In any surgical procedure the potential for adverse effects exists. Such adverse effects could include:

- Abnormal sensation and/or pain due to the implant.
- Metal sensitivity.
- Infection.
- Loosening of the wires.
- Other conditions that may place the patient at risk.

WARNINGS AND PRECAUTIONS

- The device sits proud of the skin so it should be advised to the patient that they should be conscious of this as any accidental impact may result in a lengthened recovery time.
- It is recommended that an open insertion method is used so as not to damage any sensitive or non-repairable soft tissue.
- The device should be placed so that it does not restrict the full range of motion of the hand, this particularly applies to the protruding thread.
- Inspect the sterile packaging before use, if the inner packaging is not intact, do not use the implant.
- Irrigate the fixator with saline solution during insertion to prevent heat build up in the K-wires.
- For best results, the K-wires should cross at the axis of the bone.
- Repeated adjustment of the K-wires will be to the detriment of the bone and the device and should be avoided unless absolutely necessary.
- It is advised not to apply too much compression as this can cause the fixator to bend, resulting in buckling and misalignment of the bone.
- The device is for single use only, do not attempt to re-use or re-sterilise the device.
- If the K-wires are cut too short, then it could be difficult to remove them at the end of the treatment.
- The patient should be instructed to report any unusual changes around the operated site to the physician. For example, signs of infection.

POTENTIAL COMPLICATIONS AND ADVERSE SIDE EFFECTS

- Excessive patient activity may cause the device to fail.
- Infection, pain or swelling may occur at the implant site.
- Excessive trimming of the K-wires after insertion may make removal of the device difficult.
- Loosening of the K-wires.

PACKAGING

The device is supplied sterile packed but this should only be assumed so long as the manufacturer's original packaging is intact. The device should only be received in an intact package. Damaged packages or products should not be used and should be discarded.

STERILITY

The device is supplied sterile, this is denoted on the packaging with the 'STERILE R' symbol. Sterility can only be assured provided that the original packaging is intact and that the device is within its use-by date. The device should not be used if sterility cannot be assured.

MATERIALS

The Osteo-x is manufactured from Stainless Steel to BS 7252:1, PEEK and carbon fibre.

ENVIRONMENTAL STORAGE/CONDITIONS

The Osteo-x is intended to be used and stored under normal operating environmental conditions.

SINGLE USE

The Osteo-x is intended for SINGLE USE ONLY and should under no circumstances be re-used or re-sterilised.

REPROCESSING OF SINGLE USE DEVICES

Devices labelled as SINGLE USE ONLY may not perform as intended by the manufacturer if re-used. Use of these devices causes irreversible changes to the micro and macro structure; consequently performance characteristics of the device will be sub-optimal if re-used. Re-use of a single use device may lead to:

- An increased risk of infection
- Material degradation
- Failure of the device to perform as intended
- Endotoxin reactions

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