EN

INSTRUCTIONS FOR USE

FR

MODE D'EMPLOI

DE

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ES

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IT

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SV

BRUKSANVISNING

PT

INSTRUÇÕES

TR

KULLANIM KILAVUZU

EN

INSTRUCTIONS FOR USE

MATERIAL

Skelifil HA Paste is a mixture having a nominal composition of 38% nanocrystalline hydroxyapatite in water. Hydroxyapatite has been subjected to extensive and successful clinical studies for many years. Hydroxyapatite is a safe and fully biocompatible material which is designed to act as an osteoconductive medium to support the ingrowth and fusion of adjacent viable bone when placed in an osseous environment. Skelifil HA Paste will fully resorb and be replaced by natural bone due to its nanocrystalline structure which is essentially identical to the mineral phase of human bone. Skelifil HA Paste can in many cases spare the patient the trauma of autograft harvesting and provides a sterile fully synthetic alternative to human or animal cadaver bone, and is free of the risk of disease transmission.

INDICATIONS-FOR-USE

Skelifil HA Paste is a synthetic osteoconductive bone paste that is intended to be used for filling and supporting the ingrowth of adjacent viable bone in bone defects that are not intrinsic to the stability of the bone structure. These defects may be located in the extremities, pelvis or dental, periodontal or mandibular areas.

Skelifil HA Paste may be used for filling of bone defects in the case of fractures, resection of benign tumours or cysts, and for filling autograft harvest sites. Skelifil HA Paste may also be used in defects of the alveolar ridge or wall inclusive extraction sockets, and for filling of periodontal bone defects. The paste may be gently pressed into the defect by hand and moulded to more accurately fit the defect if and as required, and may also be mixed with bone marrow aspirate prior to implantation. Skelifil HA Paste will act as a temporary scaffold and will remain as a soft paste and is not intended to provide structural support during the healing process. The implant is biocompatible, radio-opaque and will resorb in a controlled way.

CONTRAINDICATIONS

Skelifil HA Paste is not designed or sold for any other use except as indicated.

Do not use Skelifil HA Paste in the presence of any contraindication.

Skelifil HA Paste is contraindicated where the device is intended to provide structural support in the skeletal system and MUST NOT BE USED WHERE THE IMPLANTATION SITE IS UNSTABLE AND NOT RIGIDLY FIXED, OR USED TO GAIN SCREW FIXATION.

Other contraindications include:

- Existing acute or chronic infections, particularly at the implantation site
- Severe vascular or neurological disease
- Poorly vascularised implantation site
- Uncontrolled diabetes
- Severe degenerative disease
- Collagen disease
- Hypercalcaemia, abnormal calcium metabolism
- Inflammatory bone disease
- Malignant tumours
- Severely impaired renal function

PRECAUTIONS

Skelifil HA Paste is for professional use only, and is only intended for use by surgeons familiar with and skilled in the techniques of bone repair and replacement.

Skelifil HA Paste is not intended for load bearing applications. It is important to ensure that the area around the implantation site be secured mechanically with rigid fixation to provide structural support and maintain the implant in a static, load free environment. Skelifil HA Paste must not be used to gain screw fixation.

It is important to maximise the contact between existing bone and the implant to ensure proper bone regeneration.

Skelifil HA Paste should only be implanted into fresh vital bone to ensure rapid revascularisation. As with any material, care should be taken to avoid the occurrence of emboli, therefore the highly pressurised application of Skelifil HA Paste into a tightly confined defect with venous or arterial access should be avoided.

The effect of Skelifil HA Paste on patients with the following conditions is not known.

- Long term infection
- Metabolic bone disease
- Radiation bone therapy
- Cardiovascular disease

The effects of Skelifil HA Paste in paediatric patients, and the effects of mixing with other substances are unknown. However no special precautions have been identified at the time of issue.

POSSIBLE COMPLICATIONS

A successful result may not be achieved in every case.

A secondary operation to remove or replace an implant may be necessary due to surgical error, specific medical conditions or device failure. Possible adverse effects may include but are not limited to:

Displacement of the Skelifil HA paste due to a load being applied

Bone deformity or non-viable bone at the implantation site

Post operative complications including infection, haematoma, oedema, swelling, and fluid accumulation, tissue thinning, and other complications that are possible with surgery.

Allergic reaction to the product

WARNINGS

The contents of the package are double-sealed and STERILE. DO NOT USE IF OPENED, PUNCTURED OR PRODUCT OR APPLICATOR APPEARS DAMAGED OR PASTE HAS DRIED. Read the expiration date before use and DO NOT USE BEYOND THE EXPIRATION DATE.

It is recommended not to open the inner pouch until a few minutes prior to implantation.

To prevent cross-contamination the Skelifil HA Paste is for SINGLE USE ONLY. Do not attempt to re-sterilise or re-use due to the risk of cross infection.

Skelifil HA Paste should only be used in procedures where the implant can be adequately contained. Skelifil HA Paste is opaque to x-rays and may hide areas under or above the implant on a radiograph.

APPLICATION

Step 1: Open both the outer and inner pouches. Remove the cap from the applicator. Skelifil HA Paste is supplied as a paste in the applicator, ready for use.

Step 2: Implant the quantity of Skelifil HA Paste as required by the defect. The maximum dose of Skelifil HA Paste is 20g per defect. The defect should be completely filled and the paste moulded to follow the natural contours of the bone. The Skelifil HA Paste in the defect should be in direct contact with all surfaces of the defect. Where a load may be applied to the implant it should be used in combination with rigid fixation devices.

Step 3: Secure the surgical site after implanting to prevent motion and any implant displacement. When excess fluid is present at the site, cauterisation, and suction may be used to reduce bleeding.

STORAGE

Store Skelifil HA Paste at between 5°C and 25°C Optimal storage conditions: 15-25°C (59-77°F).

Direct contact with sunlight or heating systems should be avoided

SHELF LIFE AND DISPOSAL

The expiration date is printed on the labelling. Do not use Skelifil HA Paste after the expiration date.

Skelifil HA Paste is environmentally friendly. No special disposal is required. The used applicator should be disposed of as clinical waste

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